

# Hygiene Pro™ CORDLESS MICROMOTOR



*Cordless Design Maximizes Comfort*

## Lightweight and Ergonomic

The Hygiene Pro handpiece combines light weight with outstanding design. Sized and balanced specifically for today's busy dental hygiene professional and weighing in at only 3.5 ounces, the Hygiene Pro is a pleasure to use, every day.



## Smooth and Powerful

Lightweight doesn't mean you sacrifice power! The Hygiene Pro features a variable speed range (from 500 - 2,500 RPM), delivering smooth performance with all the torque you demand.

## Portable and Dependable

You will appreciate the cordless feature of the Hygiene Pro. No more drag and no more tangled cords. Plus, the Hygiene Pro's powerful battery and charger means the unit stands at the ready all day.

## Versatile

Whether you use disposable prophy angles or use screw-in or snap-on prophy cups, the Hygiene Pro has an attachment that fits your needs.

### Attachments



#### PR-F

- For Disposable Prophy Angles
- 4:1 Reduction




#### ARS-F4R

- For screw-in prophy cups
- 4:1 Reduction



#### ARK-F4R

- For snap-on prophy cups
- 4:1 Reduction



PR-F 4:1 Head for disposable prophylaxis angles

On/Off Button

Battery Life Indicator

9 Speed Settings with Indicator Light

Easy Speed Setting and Power On/Off Button