

# BiOX™

*Portable X-ray*

Innovation in X-ray imaging



*Simple Operation*

*Lightweight - 4.8lb (2.2kg)*

*Safe - Zero Leakage design  
makes BiOX safe*



**innoDEN™**  
INNOVATION MADE SIMPLE

# BiOX™

Portable X-ray

## Innovation in Today's Dental X-ray

### Portable X-ray Unit

#### • Why BiOX? •

- Lightweight! Only 4.8lbs in total weight!
- The lightest portable X-ray unit in the market today.
- ZERO Leakage design makes BiOX safe!
- One BiOX can replace all wall mounted units.
- Wide color LCD and simple operation.
- Ultra portable - You can take BiOX anywhere.



### Ergonomic Design

#### Lightweight!

The battery pack is separated from the body to minimize the actual weight of the body and to **reduce the stress on operator's arms and shoulders.**



## Recommendations

*"Using BiOX, the portable X-ray, is one of the wisest choice I've made on buying a dental product. It saves so much time that you have more consultation time with your patients and even for treatments."*

- Julie Hong, DDS

*"BiOX is simply the best portable X-ray unit, and it is the best investment for the money..."*

- Dong H. Lee, DDS

*"I prefer the BiOX over the wall mounted because it is so much lighter and easier to maneuver. It takes less than 3 minutes for me to take a FMX on a patient now and I would never want to go back to the wall mounted X-ray."*

- Jennifer Barkhouse, LDA



## Clear Image



**BiOX** produces highest quality image with lowest patient dosage. Works perfectly with both digital sensors and films resulting in excellent image. Easy and the fastest way to get the best quality x-ray image.

## Simple Operation!

**Wide color LCD** with easy to recognize preset exposure time makes **BiOX** very **easy to operate**.

You may also adjust preset exposure time by just a touch of a button.

## Affordable!

**BiOX** costs thousands of dollars less than **other competitors' portable X-ray units**. BiOX offers innovative technology and world-class quality at an affordable price. We are dedicated to make innovation simple and affordable

Ask your local dealers about BiOX. You can visit our website [www.innoden.com](http://www.innoden.com) or call 1-888-802-7379 for more information

**innoden**  
[www.innoden.com](http://www.innoden.com)



# BIOX™

Portable X-ray

## BiOX Specification

X-ray Insulation	Tube Voltage	60kV
	Tube Current	3.0mA
X-ray Tube	Anode Type	Stationary
	Anode Angle	20
	Focal Spot Size	0.8mm
Exposure Time		0.01~1.6sec (0.01 step)
Total Filtration		1.6mm Al
SSD		100mm, 210mm (W/Long Cone)
Total Weight		2.2kg (4.8lbs)
Weight W/O Battery		1.9kg (4.2lbs)
Dimension		180(W) x 73(D) x 138(H) mm
Power Requirement	Battery	CD 24V
	Battery Charger	AC 100~240V
	Power Supply (used w/stand)	AC 100~240V

